Instructions for use/
Instrucciones de uso

IMMUVIEW®

S. pneumoniae and *L. pneumophila*
Urinary Antigen Test

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EN/ES

Combined lateral flow test for qualitative detection of
S. pneumoniae and *L. pneumophila* antigens in urine
and *S. pneumoniae* in cerebrospinal fluid

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For *in vitro* diagnostic use

Intended use

The ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test is an *in vitro* rapid lateral flow test, also known as a lateral flow immunochromatographic assay. The assay is intended for qualitative detection of *Streptococcus (S.) pneumoniae* and *Legionella (L.) pneumophila* antigens in urine specimens from patients with symptoms of pneumonia. The ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test can be read visually or used in conjunction with the ImmuView® reader.

The assay is intended to aid in diagnosis of *S. pneumoniae* and of *L. pneumophila* serogroup 1 infections. The assay is further intended to aid in the diagnosis of *S. pneumoniae* infections by detection of *S. pneumoniae* antigen in cerebrospinal fluid (CSF). Results from the ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test should be interpreted in conjunction with the patient's clinical evaluation and other diagnostic methods.

Description

ImmuView *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test is a rapid lateral flow test for qualitative detection of *S. pneumoniae* in human urine and CSF samples, and *L. pneumophila* (primarily serogroup 1) antigens in human urine samples. The test is effective in presumptive diagnosis of pneumococcal pneumonia caused by *S. pneumoniae* or *Legionella* pneumonia (Legionnaires' disease) caused by *L. pneumophila*, in conjunction with culture and other

methods. Correct and early treatment is vital for the prognosis of both diseases and therefore quick methods to confirm both diseases in the initial phase are very important in order to initiate the proper treatment as soon as possible. ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test is for use by laboratory professionals and/or healthcare professionals only.

Principle

ImmuView *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test is a rapid lateral flow test for detection of *S. pneumoniae* and *L. pneumophila* using the same test.

Precautions

- Ensure that the running buffer is added to all the test tubes. False positive results can occur if no running buffer is added to the test tubes.
- The presence of partial lines and dots represents invalid test results. The patient sample should be re-tested.
- Test results should be read within the recommended reading time.
- Let the kit components equilibrate to room temperature before testing.
- Three grey/purple lines indicate plasma/protein present in the urine. The sample should be boiled for 5 minutes before retesting the sample.
- The intensity of a test line is not related to the antigen level in the sample.

- Do not mix the components of the specific kit lot with components from other kits.
- Do not use the ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test after the expiry date.
- Inspect the tests and vials before use to ensure they are intact. Any damaged vials/tests should be discarded.
- Thioglycolate as a preservative can result in a cross-reaction, potentially causing false positive *S. pneumoniae* results.

Materials provided

- 1 tube with 22 tests
- 0.5 mL combined positive control for *S. pneumoniae* and *L. pneumophila*
- 0.5 mL combined negative control for *S. pneumoniae* and *L. pneumophila*
- 2.5 mL running buffer
- 1 tweezer
- 22 transfer pipettes
- 22 test tubes
- 1 cardboard test tube holder
- 1 scorecard

Quick guide can be found on the inside of the box and on page 9.

Materials required but not provided

- Timer
- Sterile standard urine or CSF collection containers/transport tubes

Storage and stability

Please find the information on the box and labels. This product does not require any additional storage conditions. Do not freeze the product.

Preservatives

The use of boric acid or piperazine-N,N'-bis (2-ethanesulfonic acid) (PIPES) do not interfere with the ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test.

Sample collection and storage

Collect the urine sample in a sterile standard container (with or without boric acid preservative). If the sample is run within 24 hours, it can be stored at room temperature. Alternatively, the sample can be stored at 2-8 °C for 1 week or frozen at -20 °C (avoid multiple freeze/thaw cycles). Make sure that samples always reach room temperature before testing. CSF samples should be tested as soon as possible after sampling or be stored frozen until testing is possible. Follow your laboratory procedures for long term storage of CSF samples.

Quality control

The positive and negative controls provided with ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test function as the kit quality control.

The positive and negative controls should follow the same procedure as if they were urine or CSF samples. The positive control should be visible at the control test line and both the *S. pneumoniae* and *L. pneumophila* test lines. The negative control should only be visible at the control line.

Before use check the vials to ensure there is no damage and/or leak. In case of damage or leak discard the vial.

Before using a new lot of a kit, or a new shipment of the same lot, or if the test is performed by a new operator, please perform quality control testing before testing clinical samples. Follow your local or state requirements for frequency of quality control testing.

Procedure

1. Bring the patient urine or CSF sample to room temperature. Whirl sample prior to testing.*
2. Apply a test tube in the cardboard holder.
3. Fill the transfer pipette with urine or CSF and add 3 drops (120 μL) of sample to the test tube (hold the pipette vertically).
4. Add 2 drops (90 μL) of running buffer to the test tube (hold the buffer bottle vertically).
5. Whirl the test tube gently.
6. Take the container with tests, open it and take out the number of tests needed, and close it firmly afterwards.
7. Insert the test into the test tube.
8. Wait 15 minutes.
9. Lift the test out of the test tube. Read the result within 5 minutes.**
10. Discard the test after interpretation of the result.

* If the urine sample contains visible blood, please confirm a positive result by boiling^{1,2} the sample at $\geq 95^\circ\text{C}$ for 5 minutes (e.g. heating block) and retest.

** Otherwise, the test result may be inaccurate.

Quick guide

Sample addition


3 drops
(120 µL)

**Add running buffer
and whirl gently**

2 drops
(90 µL)


**Add test and wait
15 minutes**

Result interpretation

A: Control
B: *Legionella*
C: *S. pneumoniae*

*** Look closely.**
The intensity of
the lines B and C
may vary from
very clear to faint.

 1	 2	 3	 4	 5	 6
<i>Legionella</i> and <i>S. pneumoniae</i> positive	<i>Legionella</i> positive	<i>S. pneumoniae</i> positive	<i>Legionella</i> and <i>S. pneumoniae</i> positive*	Three grey/purple lines - boil and retest	Negative

Invalid test

7

No control -
retest sample

8

No control -
retest sample

9

Incomplete line -
retest sample

Interpretation of results

The control line in the top of the test will appear purple/grey but can also be more blue or red depending on whether the sample is positive for either *S. pneumoniae* or *L. pneumophila*.

A positive sample for *S. pneumoniae* and *L. pneumophila*:

A pink/red line for *S. pneumoniae* and a blue line for *L. pneumophila* and a purple/grey control line (see test results number 1 and/or number 4, page 9).

A positive sample for *L. pneumophila*:

A blue line and a purple/grey control line at the top of the test (see test result 2, page 9). A positive result for *L. pneumophila* in CSF should be investigated further, if repeatedly positive for *Legionella*.

A positive sample for *S. pneumoniae*:

A pink/red line and a purple/grey control line at the top of the test (see test result 3, page 9).

Look closely. Even if there is a very faint line for either *Legionella* or *S. pneumoniae* or both, the test result is positive (see test result number 4, page 9). The enclosed “Scorecard” can help to determine if the test result is positive or negative.

Note: Three grey/purple test lines do not indicate a positive result (see test result number 5, page 9). There can be colour differences between line A (the control line) and B/C (*S. pneumoniae*/*L. pneumophila*). If line B and C are the same colour, boil for approximately five (5) minutes and retest. Remember to let the urine sample cool down to room temperature before retesting the sample.^{1,2}

Boiling can also be used for confirmation of a positive result as *Legionella* and *S. pneumoniae* antigens are heat stable.

A negative sample:

A purple/grey control line at the top of the test (see test result 6, page 9).

Invalid sample

If no control line is observed the test is invalid and the sample should be retested (see test results number 7 and 8, page 9).

Only a full line indicates a positive result - dots and half lines do not indicate a positive result (see test result number 9, page 9).

Disposal

Follow local procedures and/or guidelines from national authorities for disposal of biological materials.

Limitations

- ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test has not been validated to be used with urine samples from children under 8 years.
- ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test has been validated using urine and CSF specimens only. Other specimens (e.g. serum, plasma or other body fluids) have not been evaluated and should not be used.

- The ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test is only validated for *S. pneumoniae* detection in CSF samples. A positive result for *L. pneumophila* in CSF should be investigated further.
- A negative result for the *Legionella* line does not exclude the possibility of a *Legionella* infection, as it can be caused by other serogroups and *Legionella* species. There is no single satisfactory laboratory test for Legionnaires' disease. Therefore, culture results, PCR, serology, and/or antigen detection methods should be used in conjunction with clinical findings to make an accurate diagnosis.
- A negative result for the *S. pneumoniae* line does not exclude an *S. pneumoniae* infection. The result of this test as well as culture, serology or other antigen detection methods should be used in conjunction with clinical findings to make an accurate diagnosis.
- *S. pneumoniae* vaccine may cause false positive results up to 10 days after vaccination.
- Administration of antibiotics might influence the test results for *S. pneumoniae*.
- False results may occur from highly basic ($\text{pH} \geq 9$) urine and give false positive *S. pneumoniae* results.
- High levels ($\geq 5\%$) of water-based personal lubricant might result in false positive or grey *L. pneumophila* lines.

Clinical sensitivity and specificity for urine samples

(Retrospective study)

The clinical sensitivity of the *S. pneumoniae* test line was obtained by testing retrospective urine samples from patients with a blood culture positive sample for *S. pneumoniae*.

The clinical sensitivity of the *L. pneumophila* test line was obtained by testing retrospective urine samples from patients with a confirmed Legionnaires' disease according to the ECDC criteria.³

The clinical specificity of the *S. pneumoniae* and *L. pneumophila* test lines was obtained by testing urine samples from patients with urinary tract infections and blood culture negative samples.³

No cross-reaction between *S. pneumoniae* and *L. pneumophila* serogroup 1 urine samples was detected. Sensitivity values were calculated using a two-sided Wilson 95% confidence interval.

ImmuView®	Sensitivity
<i>S. pneumoniae</i>	85% (60/71 CI: 74-91%)
<i>L. pneumophila</i>	89% (88/99 CI: 81-94%)
Comparator	Sensitivity
<i>S. pneumoniae</i>	77% (55/71 CI: 66-86%)
<i>L. pneumophila</i>	72% (71/99 CI: 62-80%)

ImmuView®	Specificity
<i>S. pneumoniae</i>	99% (75/76 CI: 93-100%)
<i>L. pneumophila</i>	100% (76/76 CI: 95-100%)
Comparator	Specificity
<i>S. pneumoniae</i>	100% (76/76 CI: 95-100%)
<i>L. pneumophila</i>	100% (76/76 CI: 95-100%)

Positive agreement with another urinary antigen test

S. pneumoniae positive agreement was made in a sample population containing blood culture positive samples. *L. pneumophila* positive agreement was made in a sample population containing culture and/or PCR positive samples. The positive agreement was calculated in accordance with a two-sided Wilson confidence interval 95%.

<i>S. pneumoniae</i> blood culture positive samples		<u>Comparator</u>		
		Positive	Negative	Total
ImmuView®	Positive	58	5	63
	Negative	2	11	13
Total		60	16	76
Positive agreement		97% (CL: 89-99%)		

<i>L. pneumophila</i> culture and/or PCR positive		<u>Comparator</u>		
		Positive	Negative	Total
ImmuView®	Positive	64	10	74
	Negative	1	8	9
Total		65	18	83
Positive agreement		98% (CL: 92-100%)		

Negative agreement with another urinary antigen test

S. pneumoniae negative agreement was made in a sample population containing non-pneumococcal cases (n = 90 bacteremic samples and n = 6 non-bacteremic). *L. pneumophila* negative agreement was made in a sample population pre-screened with urinary antigen EIA. All samples included in the pool showed a negative result when using the comparator (urinary antigen EIA).

96 non-pneumococcal cases		Comparator		
		Positive	Negative	Total
ImmuView®	Positive	3	0	3
	Negative	0	93	93
Total		3	93	96
Negative agreement		100% (CL: 96-100%)		

Negative urinary antigen EIA directed at <i>L. pneumophila</i> sg. 1		Comparator		
		Positive	Negative	Total
ImmuView®	Positive	0	0	0
	Negative	0	456	456
Total		0	456	456
Negative agreement		100% (CL: 96-100%)		

Analytical studies - urine

Specificity (cross-reactivity)

To determine the analytical specificity of the ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test urines were spiked with whole cell bacteria and different inactivated viruses (N=143). The whole cell bacterial panel was tested in a concentration of 10^7 CFU/mL diluted from a stock solution. The viral panel had a concentration of 10^5 TCID₅₀/mL. The panel was also tested in negative urine.

Organisms tested for interference

<i>Acinetobacter</i> spp. (4)	<i>Lactobacillus</i> sp.
<i>Bacillus subtilis</i>	<i>Listeria monocytogenes</i>
<i>Bordetella pertussis</i>	<i>Morganella morganii</i>
<i>Moraxella catarrhalis</i>	<i>Moraxella osloensis</i>
<i>Candida albicans</i> (4)	<i>Mycoplasma genitalium</i>
<i>Citrobacter freundii</i>	<i>Neisseria gonorrhoeae</i> (3)
<i>Cornyebacterium</i> sp.	<i>Neisseria lactamica</i>
<i>Cornyebacterium uralyticum</i>	<i>Neisseria meningitidis</i>
<i>Enterobacter cloacae</i> (3)	<i>Neisseria polysaccharea</i>
<i>Escherichia coli</i> (10)	<i>Proteus mirabilis</i> (2)
<i>Enterococcus faecalis</i> (7)	<i>Proteus vulgaris</i>
<i>Enterococcus faecium</i>	<i>Pseudomonas aeruginosa</i> (4)
<i>Enterococcus durans</i>	<i>Pseudomonas stutzeri</i>
<i>Gardnerella vaginalis</i>	<i>Pseudomonas</i> spp. (2)
<i>Haemophilus influenzae</i> type a-f and non-caps (11)	<i>Salmonella</i> Bredeney
<i>Haemophilus parainfluenzae</i>	<i>Salmonella</i> Thompson
Adenovirus 2	<i>Salmonella</i> Typhimurium
<i>Chlamydophila pneumoniae</i> (2)	<i>Serratia marcescens</i>
<i>Chlamydia trachomatis</i>	<i>Staphylococcus epidermidis</i>
Cytomegalovirus	<i>Salmonella</i> Glostrup
Enterovirus D68	<i>Streptococcus mutans</i> (2)
Herpes Simplex 1,2	<i>Streptococcus parasanguis</i>

Organisms tested for interference	
Influenza A (H1N1 and H3N2) virus	<i>Streptococcus sanguinis</i>
Influenza B Virus	<i>Streptococcus aureus</i> (6)
Parainfluenza virus 1,2,3 (3)	<i>Streptococcus epidermidis</i> (5)
Respiratory Syncytial Virus A	<i>Streptococcus saprophyticus</i> (3)
<i>Klebsiella oxytoca</i> (2)	<i>Stenotrophomonas maltophilia</i>
<i>Klebsiella pneumoniae</i> (3)	<i>Streptococcus</i> Gr. A, B, C, F, L and G (16)
<i>Lactobacillus cateniforme</i>	<i>Streptococcus mitis</i>
<i>Lactobacillus rhamnosus</i>	

All of the above bacterial isolates were negative when using ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test. The only potential cross-reactivity was 1 of 3 isolates of *E. cloacea* which was positive for *L. pneumophila*. This was confirmed on re-testing of that one isolate.

A total of 19 urine samples from patients with urinary tract infection were tested. Previously, culture results had shown that eight (8) of them were infected with *Escherichia coli*, five (5) with *Staphylococcus aureus*, five (5) with *Streptococcus agalactiae* Gr. B and one (1) with *Candida albicans*. None showed any cross reactions with the ImmuView® test.

Sensitivity (limit of detection (LoD))

The limit of detection for the ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test is 62.5 pg/mL for purified *S. pneumoniae* CWPS antigen (native). For LPS specific for *L. pneumophila* SG1 (Philadelphia) the limit of detection is 25 ng/mL. Whole cell *S. pneumoniae* bacteria can be detected at a limit of detection at 10^5 CFU/mL and *L. pneumophila* SG1 (Philadelphia) has a limit of detection at 10^4 CFU/mL. Boiling or urine preservatives did not change these results.

Content	LoD
<i>S. pneumoniae</i> antigen	62.5 pg/mL
<i>L. pneumophila</i> SG 1 (Philadelphia) antigen	0.025 µg/mL
<i>L. pneumophila</i> SG 1 (Bellingham) antigen	0.5 µg/mL
<i>S. pneumoniae</i> (serotype 1)	10^5 CFU/mL
<i>L. pneumophila</i> SG1 (Philadelphia)	10^4 CFU/mL
<i>L. pneumophila</i> SG 1 (Bellingham)	10^5 CFU/mL

Strain reactivity

Isolates from different *S. pneumoniae* serotypes were also tested positive with the ImmuView® assay including serotype three (3), five (5) and thirty-seven (37). Different species of *L. pneumophila* were also found to be positive using the assay. Within serogroup one (1) these includes Philadelphia, Knoxville, OLDA/Oxford, Allentown/France, and Benidorm-Strain Lens. In addition, previous studies have found other *Legionella* serogroups to be positive such as serogroup three (3), six (6) and twelve (12).⁴

<i>S. pneumoniae</i>		
Serotype	Antigen (µg/mL)	Whole organism (CFU/mL)
1	ND	10 ⁴
3	0.001	10 ⁴
5	0.010	10 ⁵
37	0.0001	ND

<i>L. pneumophila</i> serogroup one (1)			
Pontiac/non-Pontiac	Species	Antigen (µg/mL)	Whole cells (CFU/mL)
Pontiac	Knoxville	0.100	10 ⁵
Pontiac	Allentown/France	0.005	ND
Pontiac	Benidorm	ND	10 ⁴
Pontiac	Philadelphia	0.010	10 ⁴
Non-Pontiac	OLDA/Oxford	0.001	ND
Non-Pontiac	Camperdown	0.315	ND
Non-Pontiac	Heyman	1.250	ND

<i>L. pneumophila</i> non-serogroup one (1) detected by ImmuView ⁴
Serogroup; 2, 3, 4, 5, 6, 7, 11, 12, 13, 14 and 15

Interfering substances

ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test was tested with forty-seven (47) interfering agents at different concentrations in urine samples.

Agent	Concentration	Agent	Concentration
Acetaminophen	0.1mg/mL	Leucocytes	>250 cells/μL
Acetylsalicylic acid	0.1mg/mL	Miconazole	5%
Amantadine	0.03mg/mL	Mix (pH, whole blood, protein and glucose)	N/A
Amoxicillin	0.075mg/mL	Mucin	0.086mg/mL
Amphotericin B	0.22mg/mL	Oseltamivir (Tamiflu)	0.03mg/mL
Antihistamine	0.22mg/mL	Oxalic acid	0.01%
Ascorbic acid (C-Vitamin)	1mg/mL	pH (acidic)	4
Augmentin (Amoxicillin Clavulanate)	0.22mg/mL	pH (neutral)	7
Azithromycin	0.012mg/mL	pH (basic)	9
Beet root (H)	20%	Plasma	90%
Beet root (M)	1.17%	Plasma	50%
Beet root (L)	0.01%	Plasma	10%
Bilirubin	0.2mg/mL	Prednisone	0.22mg/mL
Bromhexin/cough drops/ cough syrup	0.22mg/mL	Protein (albumin) (H)	10mg/mL
Caffeine	15mg/mL	Protein (albumin) (M)	5mg/mL
Chlorophyll (H)	0.11mg/mL	Protein (albumin) (L)	0.6mg/mL
Chlorophyll (M)	0.04mg/mL	Pyridium	1mg/mL
Chlorophyll (L)	0.01mg/mL	Rifampicin	0.09mg/mL
Ciprofloxacin	0.22mg/mL	Spinach	1%
Decongestant	0.22mg/mL	Tobacco purified	0.4mg/mL
Corticosterone (Corticosteroids)	0.015mg/mL	Triglycerides	4mg/mL
Erythromycin	0.067mg/mL	Urea	20mg/mL
Glucose (H)	20mg/mL	Vaginal contraceptive gel	5%
Glucose (M)	10mg/mL	Vancomycin	0.1mg/mL
Glucose (L)	3mg/mL	Water-based personal lubricant	5%
Hemoglobin	5mg/mL	White blood cells	10%
Human albumin	35mg/mL	Whole blood	10%
Human red blood cells 10% Washed pooled cells	10%	Whole blood	15%
Ibuprofen	0.1mg/mL		
Itraconazole	0.22mg/mL		

High concentration of plasma in urine may result in grey test lines. Additionally, basic ($\text{pH} \geq 9$) conditions in urine can give false positive *S. pneumoniae* lines. Water-based personal lubricant might result in false positive or grey *L. pneumophila* lines, however, this outcome seems dose-related.

Clinical sensitivity and specificity for CSF

The sensitivity of the *S. pneumoniae* was obtained by testing 12 CSF samples which were culture positive, and 15 CSF samples spiked with *S. pneumoniae*. The specificity of the *S. pneumoniae* and *L. pneumophila* test lines were obtained by testing 170 negative CSF samples from negative donors.

ImmuView®	Sensitivity	Specificity
<i>S. pneumoniae</i>	100% (27/27)	98.8% (168*/170)
<i>L. pneumophila</i>	N/A	100% (170/170)

* 2 samples were tested positive and confirmed positive with both another lateral flow test for *S. pneumoniae* and ImmuLex *S. pneumoniae* Omni. It was not possible to culture any bacteria from the samples, which can be caused by too many times of freezing and thawing of the sample.

Spiked CSF testing

50 human CSF samples spiked to the limit of detection with *S. pneumoniae* and 10 negative CSF samples were tested with the ImmuView® test and the comparator test. The sensitivity for the ImmuView® test and the comparator test was 50/50 (100%) whereas the negative CSF samples showed 10/10 (100%) for both tests.

60 real human CSF samples			
ImmuView®	<u>Comparator</u>		Total
	Positive	Negative	
Positive	50	0	50
Negative	0	10	10
Positive agreement	100%	95% CI (92.9%-100%)	
Negative agreement	100%	95% CI (72.2%-100%)	

Analytical studies - CSF

Specificity (cross-reactivity)

ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test was tested with a panel of 24 potential cross-reacting agents. No cross-reactions were detected for the *S. pneumoniae* or the *L. pneumophila* test lines.

Organisms not affecting test performance in CSF	
<i>Escherichia coli</i> (5)	<i>Neisseria meningitidis</i> Gr. B, D and W135 (3)
<i>Haemophilus influenzae</i> type a-f and non-caps (7)	<i>Staphylococcus aureus</i>
<i>Listeria monocytogenes</i>	<i>Streptococcus</i> Gr A
Measles	<i>Streptococcus agalactiae</i> (GBS) sg Ia, Ib, II, III (4)
<i>Streptococcus mitis</i>	

Sensitivity (limit of detection (LoD)) in CSF

The analytical sensitivity for the ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test was determined by limit of detection. Two different operators performed the dilutions and the testing. The dilutions were made with whole cell bacteria spiked in human CSF.

CSF	LoD
<i>S. pneumoniae</i>	10 ³ CFU/mL

Interfering substances

ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test was tested with forty-seven (47) interfering agents at different concentrations in artificial CSF either negative or spiked with either CWPS (10ng/mL) or *S. pneumoniae* 10⁷ CFU/mL.

Agent in CSF	Concentration	Agent	Concentration
<i>S. pneumoniae</i> (Type 1) 10⁷ CFU/mL		Negative artificial CSF	
Glucose (H)	1mg/mL	Glucose (H)	1mg/mL
Glucose (M)	0.5mg/mL	Glucose (M)	0.5mg/mL
Glucose (L)	0.1mg/mL	Glucose (L)	0.1mg/mL
Red blood cells (H)	15%	Red blood cells (H)	15%
Red blood cells(M)	10%	Red blood cells(M)	10%
Red blood cells (L)	5%	Red blood cells (L)	5%
Protein (H)	60mg/mL	Protein (H)	60mg/mL
Protein (M)	30mg/mL	Protein (M)	30mg/mL
Protein (L)	10mg/mL	Protein (L)	10mg/mL
White blood cells	10.6x10 ⁶ /mL	White blood cells	10.6x10 ⁶ /mL
White blood cells	5.3x10 ⁶ /mL	White blood cells	5.3x10 ⁶ /mL
White blood cells	2.7x10 ⁶ /mL	White blood cells	2.7x10 ⁶ /mL
White blood cells	1.8x10 ⁶ /mL	White blood cells	1.8x10 ⁶ /mL
White blood cells	0.9x10 ⁶ /mL	White blood cells	0.9x10 ⁶ /mL
Antigen (CWPS) 10 ng/mL		Bilirubin	15%
Bilirubin	15%	Bilirubin	10%
Bilirubin	10%	Bilirubin	5%
Bilirubin	5%	Plasma	15%
Plasma	15%	Plasma	10%
Plasma	10%	Plasma	5%
Plasma	5%		

Red blood cells may give false positive shadows on the *S. pneumoniae* line due to excessive red colour. The other agents in the panel did not interfere with the test.

Reproducibility study

The ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen test demonstrated excellent overall reproducibility with 1,068 correct results out of 1,072 test results (99.6%), when tested with 10 clinically positive *S. pneumoniae* or *L. pneumophila* urine samples and negative urine samples, and artificial CSF positive spiked with *S. pneumoniae* isolates as well as negative artificial CSF samples. The ImmuView® positive control and negative control were also tested as blinded/masked panel members. The testing was performed for 5 days with a different kit lot at each site, two in the U.S. and one in Europe.

A total of 3 different lots were tested at the 3 different sites using two operators each (A and B). A total of 8 test results (0.7%) were determined to be invalid and were excluded and not re-tested. The panel members were blinded by changing the panel member numbers and identity daily. The reading and interpretation of the reproducibility panels were performed visually. There were no statistical differences in reproducibility by lot, by site, by time or by operator.

Description	Correct results	Agreement
<i>S. pneumoniae</i> , moderate positive urine	90/90 positive	100.0%
<i>S. pneumoniae</i> , moderate positive CSF	89/89 ^A positive	100.0%
<i>S. pneumoniae</i> , low positive spiked in artificial CSF	89/90 ^B positive	98.9%
<i>S. pneumoniae</i> , low positive urine	90/90 positive	100.0%
<i>L. pneumophila</i> , moderate positive urine 2A	90/90 positive	100.0%
<i>L. pneumophila</i> , moderate positive urine 2B	88/89 ^C positive	98.9%
<i>L. pneumophila</i> , low positive urine 1A	89/89 ^D positive	100.0%
<i>L. pneumophila</i> , low positive urine 1B	89/90 ^E positive	98.9%
Negative pooled urine	90/90 negative	100.0%
Negative artificial CSF	90/90 negative	100.0%
ImmuView [®] positive control	89/90 ^F positive	98.9%
ImmuView [®] negative control	85/85 ^G negative	100.0%
Summary	1068/1072 correct	99.6%

- A Operator did not see a positive control band, so one sample was invalid as the package insert states that this is necessary before interpreting the result. The sample was not re-tested.
- B A visual *L. pneumophila* band was seen.
- C Operator interpreted band as *S. pneumoniae* positive instead of *L. pneumophila* positive. One sample was invalid due to dots/incomplete band on the test per the package insert and was not re-tested.
- D One sample was invalid due to an incomplete band in *S. pneumoniae* according to the pack insert.
- E No *L. pneumophila* band present.
- F Operator interpreted *S. pneumoniae* band result as negative even though band was present.
- G Five samples excluded due to the presence of dots/incomplete bands. The samples were not re-tested.

Incident reporting

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the member state in which the user and/or patient is established.

Quality certificate

SSI Diagnostica's development, production and sales of *in vitro* diagnostics are quality assured and certified in accordance with ISO 13485.

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SSI Diagnostica A/S