

## General data sheet for transport and storage media

### Intended use

Transport and storage media are used to preserve microorganisms alive and enable transport and possible propagation of microorganisms.

For use by clinical laboratories as part of their established procedures.  
Analytical and clinical performance have not been determined by the manufacturer.

### Principle

The contents of the various tubes are adapted to the needs of the individual microorganisms. They can be used for propagation of the microorganism or for storage.

### Quality control

The quality control may consist of chemical and or physical testing, depending on the use of the product.

### Quality control at user

The user is advised to perform quality control depending on the use of the product.

### Storage and stability

See storage and shelf life on the label.

### Precautions

For professional use only for *in vitro* diagnostics. IVD

Do not use the product if there are signs of microbial contamination, discolouration, drying out or other signs of obsolescence. Sown products must be disposed of in accordance with current regulations on clinical risk waste.

### Quality assurance

SSI Diagnostica is quality assured and certified in accordance with ISO 13485. Analysis certificates and safety data sheets can be downloaded from our website [www.ssidiagnostica.com](http://www.ssidiagnostica.com).

**Information and ordering**

SSI Diagnostica A/S

Herredsvejen 2

3400 Hillerød

Denmark

T +45 4829 9100

[info@ssidiagnostica.com](mailto:info@ssidiagnostica.com)

[ssidiagnostica.com](http://ssidiagnostica.com)

[shop.ssidiagnostica.com](http://shop.ssidiagnostica.com)