


FLEXICULT® URINKIT

Escherichia coli

ATCC 35218




10³ CFU/mL



10⁴ CFU/mL




10⁵ CFU/mL




Klebsiella pneumoniae

ATCC 13883





10³ CFU/mL




10⁵ CFU/mL

Enterobacter cloacae

ATCC 13047



10³ CFU/mL



10⁵ CFU/mL

Proteus mirabilis

ATCC 12453



10³ CFU/mL



10⁵ CFU/mL

Proteus vulgaris

ATCC 8427



10³ CFU/mL



10⁵ CFU/mL

Pseudomonas aeruginosa

ATCC 27853




10³ CFU/mL



10⁵ CFU/mL

Enterococcus faecalis

ATCC 29212



10³ CFU/mL



10⁵ CFU/mL

Enterococcus faecium

ATCC 35667



10³ CFU/mL



10⁵ CFU/mL

Staphylococcus saprophyticus

ATCC 15305



10³ CFU/mL



10⁵ CFU/mL

Streptococcus agalactiae

ATCC 12927



10³ CFU/mL



10⁵ CFU/mL

Candida albicans



ATCC 2091



10⁵ CFU/mL

Mixed Flora *Escherichia coli* / *Enterococcus faecalis*

ATCC 35218 / ATCC 29212




10³ CFU/mL



10⁵ CFU/mL

Flexicult® HUMAN
contains the following
antibiotics:
1 = Mecillinam
2 = Nitrofurantoin
3 = Ampicillin
4 = Sulfamethizol
5 = Trimethoprim



All the agar plates are photographed
from above without the lid

For more information please
see the Instruction for Use.